



**Acuity Size**

0 20/400  
 9 20/200  
 1 20/100  
 8 20/80  
 7 20/70  
 6 20/60  
 5 20/50  
 4 20/40  
 3 20/30  
 2 20/20  
 . Reset

**Controls**

H Hide/unhide  
 L Single line  
 U Multi line  
 V Video  
 Space Randomize  
 R, \* Temp. Random  
 F4 Random on/off  
 F1 Help Key  
 O One Letter

**Acuity type**

S Snellen  
 I ETDRs  
 P Pictures  
 E Tumbling E  
 C Tumbling C  
 N Numbers  
 T HOTV

**Special test**

W Worth dot  
 F2 Fixations disparity  
 A Astigmatic Display  
 J Reduced Contrast  
 G Red/Green  
 K Crowding Bar

Page Down, Decrease  
 - acuity 1  
 line  
 Page up Increase  
 - acuity 1  
 line